SPORT-ELEC®

I N S T I T U T

FR

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GB

UTILITY OF PRODUCT

- Treatment of traumas
- Inflammation
- Arthritis
- Post-op protections
- Rehabilitative applications, in case it needs compression and warmth to the knee, with ligament and joint problems, and also efficient for osteoarthritis.

MATERIAL

- Neopren group: 5 mm (neoprene 80%, cotton 15%, polyester 5%).
- Woven group: Lycra woven textile (55% polyamide, 45% elastan).

IMPORTANT ADVISES AND WARNINGS

- This product is not burn protected. It will lose its characteristics by high temperature. Therefore do not use above 248°F.
- If you want to make changes for your comfort, please ask to your doctor or where you purchased the product.
- This product can not be used out of its own utility and cannot be combined with other products, only use with its original utilities.
- Usage of neoprene products is not allowed for more than 3-4 hours.
- In any case, do not use more than 8 hours (unless otherwise specified by your doctor) Please only use with doctor advice.

WASHING INSTRUCTIONS

- Please close the velcro bands if there are and take out the locable support (like joint parts...).
- Do not put in washing-machine. Only handwash at 104°F.
- Do not use bleach.
- No dry-cleaning.
- Do not use dryer.
- No ironing.
- After the wash please be sure that it is totally dry before the use.

UTILITY OF THE PRODUCT

- Please be sure the size of the product is exactly your size.
- Please open the velcro band. Be sure the kneepad completely wraps the patella. Feel that the tip of the patella is free.

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